

IL13RA2 antibody

Product Information

Catalog No.:	FNab04221
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Binds as a monomer with high affinity to interleukin-13(IL13), but not to interleukin-4(IL4).

Immunogen information


Immunogen:	interleukin 13 receptor, alpha 2
Synonyms:	Interleukin-13 receptor subunit alpha-2 (IL-13 receptor subunit alpha-2, IL-13R subunit alpha-2, IL-13R-alpha-2, IL-13RA2) Interleukin-13-binding protein IL13RA2 IL13R
Observed MW:	50-65 kDa
Uniprot ID :	Q14627

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab04221(IL13RA2,CD213A2 Antibody) at dilution of 1:200



rat thymus tissue were subjected to SDS PAGE followed by western blot with FNab04221(IL13RA2,CD213A2 Antibody) at dilution of 1:300