

IL1F10 antibody

Product Information

Catalog No.:	FNab04236
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	9H1
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

Cytokine with immunomodulatory activity. Alone, does not induce cytokine production, but reduces IL22 and IL17A production by T-cells in response to heat-killed *Candida albicans*. Reduces IL36G-induced production of IL8 by peripheral blood mononuclear cells. Increases IL6 production by dendritic cells stimulated by bacterial lipopolysaccharides(LPS). Ligand for IL-36R/IL1RL2.


Immunogen information

Immunogen:	interleukin 1 family, member 10(theta)
Synonyms:	Interleukin-1 family member 10 (IL-1F10) Family of interleukin 1-theta (FIL1 theta) Interleukin-1 HY2 (IL-1HY2) Interleukin-1 theta (IL-1 theta) Interleukin-38 (IL-38) IL1F10 FIL1T IL1HY2 IL38
Observed MW:	17 kDa
Uniprot ID :	Q8WWZ1

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000

Image:



Recombinant protein were subjected to SDS PAGE followed by western blot with FNab04236(IL1F10 Antibody) at dilution of 1:1000