

IL21R antibody

Product Information

Catalog No.:	FNab04252
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


This is a receptor for interleukin-21.

Immunogen information


Immunogen:	interleukin 21 receptor
Synonyms:	Interleukin-21 receptor (IL-21 receptor, IL-21R) Novel interleukin receptor IL21R NILR
Observed MW:	68-70 kDa
Uniprot ID :	Q9HBE5

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lymphoma using FNab04252(IL21R antibody) at dilution of 1:10



Raji cells were subjected to SDS PAGE followed by western blot with FNab04252(IL21R antibody) at dilution of 1:500