

## IL24 antibody

### Product Information


Catalog No.:	FNab04257
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	interleukin 24
Synonyms:	Interleukin-24 (IL-24) Melanoma differentiation-associated gene 7 protein (MDA-7) Suppression of tumorigenicity 16 protein IL24 MDA7 ST16
Observed MW:	24 kDa
Uniprot ID :	Q13007

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse lung were subjected to SDS PAGE followed by western blot with FNab04257(IL24 Antibody) at dilution of 1:600