

IL27 antibody

Product Information

Catalog No.:	FNab04260
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	1H3
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	interleukin 27
Synonyms:	Interleukin-27 subunit alpha (IL-27 subunit alpha, IL-27-A, IL27-A) Interleukin-30 p28 IL27 IL27A IL30
Observed MW:	31 kDa
Uniprot ID :	Q8NEV9

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:2000-1:10000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta tissue slide using FNab04260(IL27 Antibody) at dilution of 1:50



IL27 fusion protein were subjected to SDS PAGE followed by western blot with FNab04260(IL27 Antibody) at dilution of 1:5000