

IL31 antibody

Product Information

Catalog No.:	FNab04269
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Interleukin-31(IL-31) is a newly discovered four-helix bundle cytokine and expressed by several types of cells under both normal and disease state, such as activated CD4+T cells, mast cells, monocytes, macrophages, immature and mature monocyte-derived dendritic cells and so on. IL-31 signals through a heterodimeric receptor complex consisting of the IL31 receptor alpha(IL31Rα) and the oncostatin

Immunogen information

Immunogen:	interleukin 31
Synonyms:	Interleukin-31 (IL-31) IL31
Observed MW:	None
Uniprot ID :	Q6EBC2

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using FNab04269(IL31 Antibody) at dilution of 1:50