

IL3RA antibody

Product Information

Catalog No.:	FNab04277
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


This is a receptor for interleukin-3.

Immunogen information

Immunogen:	interleukin 3 receptor, alpha(low affinity)
Synonyms:	Interleukin-3 receptor subunit alpha (IL-3 receptor subunit alpha, IL-3R subunit alpha, IL-3R-alpha, IL-3RA) IL3RA IL3R
Observed MW:	70 kDa
Uniprot ID :	P26951

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HL-60 cells were subjected to SDS PAGE followed by western blot with FNab04277(IL3RA antibody) at dilution of 1:300