

## IMP4 antibody

### Product Information

Catalog No.:	FNab04298
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Component of the 60-80S U3 small nucleolar ribonucleoprotein(U3 snoRNP). Required for the early cleavages during pre-18S ribosomal RNA processing.

### Immunogen information

Immunogen:	IMP4, U3 small nucleolar ribonucleoprotein, homolog(yeast)
Synonyms:	U3 small nucleolar ribonucleoprotein protein IMP4 (U3 snoRNP protein IMP4) Brix domain-containing protein 4 IMP4 BXDC4
Observed MW:	34 kDa
Uniprot ID :	Q96G21

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04298(IMP4 antibody) at dilution of 1:500