

INPP5B antibody

Product Information

Catalog No.:	FNab04327
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Hydrolyzes phosphatidylinositol 4,5-bisphosphate(PtIns(4,5)P₂) and the signaling molecule phosphatidylinositol 1,4,5-trisphosphate(PtIns(1,4,5)P₃), and thereby modulates cellular signaling events.

Immunogen information


Immunogen:	inositol polyphosphate-5-phosphatase, 75kDa
Synonyms:	Type II inositol 1,4,5-trisphosphate 5-phosphatase 75 kDa inositol polyphosphate-5-phosphatase Phosphoinositide 5-phosphatase (5PTase) INPP5B OCRL2
Observed MW:	75 kDa
Uniprot ID :	P32019

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lymphoma using FNab04327(INPP5B antibody) at dilution of 1:50



mouse spleen tissue were subjected to SDS PAGE followed by western blot with FNab04327(INPP5B antibody) at dilution of 1:200