

ITGA3 antibody

Product Information

Catalog No.:	FNab04345
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	4F6
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


ITGA3, also named as MSK18, CD49c, GAPB3 and FRP-2, belongs to the integrin alpha chain family. It can be cleaved into heavy chain(840aa) and light chain(176aa). This antibody can recognize the heavy chain.

Immunogen information

Immunogen:	integrin, alpha 3
Synonyms:	Integrin alpha-3 CD49 antigen-like family member C FRP-2 Galactoprotein B3 (GAPB3) VLA-3 subunit alpha Integrin alpha-3 heavy chain Integrin alpha-3 light chain ITGA3 MSK18
Observed MW:	120 kDa
Uniprot ID :	P26006

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab04345(ITGA3 antibody) at dilution of 1:1000