

## INTS12 antibody

### Product Information

Catalog No.:	FNab04360
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Component of the Integrator complex, a complex involved in the small nuclear RNAs(snRNA) U1 and U2 transcription and in their 3'-box-dependent processing. The Integrator complex is associated with the C-terminal domain(CTD) of RNA polymerase II largest subunit(POLR2A) and is recruited to the U1 and U2 snRNAs genes.

### Immunogen information

Immunogen:	integrator complex subunit 12
Synonyms:	Integrator complex subunit 12 (Int12) PHD finger protein 22 INTS12 PHF22
Observed MW:	49-52 kDa
Uniprot ID :	Q96CB8

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab04360(INTS12 antibody) at dilution of 1:1000