

IQGAP2 antibody

Product Information

Catalog No.:	FNab04379
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Binds to activated CDC42 and RAC1 but does not seem to stimulate their GTPase activity.
Associates with calmodulin.

Immunogen information


Immunogen:	IQ motif containing GTPase activating protein 2
Synonyms:	Ras GTPase-activating-like protein IQGAP2 IQGAP2
Observed MW:	181 kDa, 65 kDa
Uniprot ID :	Q13576

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab04379(IQGAP2 Antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab04379(IQGAP2 Antibody) at dilution of 1:300