

## IQGAP3 antibody

### Product Information

Catalog No.:	FNab04380
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	IQ motif containing GTPase activating protein 3
Synonyms:	Ras GTPase-activating-like protein IQGAP3 IQGAP3
Observed MW:	185 kDa
Uniprot ID :	Q86VI3

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab04380(IQGAP3 antibody) at dilution of 1:2000