

IRF5 antibody

Product Information

Catalog No.:	FNab04392
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	7C4
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Transcription factor involved in the induction of interferons IFNA and INFB and inflammatory cytokines upon virus infection. Activated by TLR7 or TLR8 signaling.

Immunogen information

Immunogen:	interferon regulatory factor 5
Synonyms:	Interferon regulatory factor 5 (IRF-5) IRF5
Observed MW:	56 kDa
Uniprot ID :	Q13568

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500 - 1:3000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04392(IRF5 antibody) at dilution of 1:1000