

ISCA2 antibody

Product Information

Catalog No.:	FNab04401
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in the maturation of mitochondrial 4Fe-4S proteins functioning late in the iron-sulfur cluster assembly pathway. May be involved in the binding of an intermediate of Fe/S cluster assembly.

Immunogen information

Immunogen:	iron-sulfur cluster assembly 2 homolog
Synonyms:	Iron-sulfur cluster assembly 2 homolog, mitochondrial HESB-like domain-containing protein 1 ISCA2 HBLD1
Observed MW:	None
Uniprot ID :	Q86U28

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



A431 cells were subjected to SDS PAGE followed by western blot with FNab04401 (ISCA2 antibody) at dilution of 1:300