

ISG20L2 antibody

Product Information

Catalog No.:	FNab04405
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


3'-> 5'-exoribonuclease involved in ribosome biogenesis in the processing of the 12S pre-rRNA. Displays a strong specificity for a 3'-end containing a free hydroxyl group.

Immunogen information

Immunogen:	interferon stimulated exonuclease gene 20kDa-like 2
Synonyms:	Interferon-stimulated 20 kDa exonuclease-like 2 ISG20L2
Observed MW:	40 kDa
Uniprot ID :	Q9H9L3

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab04405(ISG20L2 Antibody) at dilution of 1:1000