

REG3G antibody

Product Information

Catalog No.:	FNab04406
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	ISL LIM homeobox 1
Synonyms:	Regenerating islet-derived protein 3-gamma (REG-3-gamma) Pancreatitis-associated protein 1B (PAP-1B) Pancreatitis-associated protein 1B (PAP 1B) Regenerating islet-derived protein III-gamma (REG III, Reg III-gamma) Regenerating islet-derived protein 3-gamma 16.5 kDa form Regenerating islet-derived protein 3-gamma 15 kDa form REG3G PAP1B
Observed MW:	39 kDa, 45 kDa
Uniprot ID :	Q6UW15

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP, IHC
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas tissue slide using FNab04406(Islet 1 Antibody) at dilution of 1:200



IP Result of anti-Islet 1 (IP:FNab04406, 4ug; Detection:FNab04406 1:500) with HeLa cells lysate 2000ug.



BxPC-3 cells were subjected to SDS PAGE followed by western blot with FNab04406(ISL1 antibody) at dilution of 1:500