

ITFG1 antibody

Product Information

Catalog No.:	FNab04410
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Modulator of T-cell function. Has a protective effect in graft versus host disease model(By similarity).

Immunogen information

Immunogen:	integrin alpha FG-GAP repeat containing 1
Synonyms:	T-cell immunomodulatory protein (Protein TIP) Integrin-alpha FG-GAP repeat-containing protein 1 Linkin ITFG1 LNKN-1 TIP
Observed MW:	75 kDa
Uniprot ID :	Q8TB96

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000
Image:	



A431 cells were subjected to SDS PAGE followed by western blot with FNab04410(ITFG1 antibody) at dilution of 1:1000