

JAKMIP2 antibody

Product Information

Catalog No.:	FNab04434
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	janus kinase and microtubule interacting protein 2
Synonyms:	Janus kinase and microtubule-interacting protein 2 CTCL tumor antigen HD-CL-04 Neuroendocrine long coiled-coil protein 1 JAKMIP2 JAMIP2 KIAA0555 NECC1
Observed MW:	95 kDa
Uniprot ID :	Q96AA8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab04434(JAKMIP2 antibody) at dilution of 1:1000