

MAPK8IP3 antibody

Product Information

Catalog No.:	FNab04438
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	mitogen-activated protein kinase 8 interacting protein 3
Synonyms:	C-Jun-amino-terminal kinase-interacting protein 3 (JIP-3, JNK-interacting protein 3) JNK MAP kinase scaffold protein 3 Mitogen-activated protein kinase 8-interacting protein 3 MAPK8IP3 JIP3 KIAA1066
Observed MW:	147 kDa
Uniprot ID :	Q9UPT6

Application

Reactivity:	Human, Rat, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



rat brain tissue were subjected to SDS PAGE followed by western blot with FNab04438(JIP3 Antibody) at dilution of 1:600