

JMJD4 antibody

Product Information

Catalog No.:	FNab04441
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	jumonji domain containing 4
Synonyms:	2-oxoglutarate and iron-dependent oxygenase JMJD4 JmjC domain-containing protein 4 Jumonji domain-containing protein 4 Lysyl-hydroxylase JMJD4 JMJD4
Observed MW:	48 kDa
Uniprot ID :	Q9H9V9

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04441(JMJD4 antibody) at dilution of 1:600