

MAPK8 antibody

Product Information

Catalog No.:	FNab04445
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	mitogen-activated protein kinase 8
Synonyms:	Mitogen-activated protein kinase 8 (MAP kinase 8, MAPK 8) JNK-46 Stress-activated protein kinase 1c (SAPK1c) Stress-activated protein kinase JNK1 c-Jun N-terminal kinase 1 MAPK8 JNK1 PRKM8 SAPK1 SAPK1C
Observed MW:	44-48 kDa, 50-55 kDa
Uniprot ID :	P45983

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human cervical cancer tissue slide using FNab04445(JNK Antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab04445(JNK antibody) at dilution of 1:1000