

MAPK9 antibody

Product Information

Catalog No.:	FNab04448
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	mitogen-activated protein kinase 9
Synonyms:	Mitogen-activated protein kinase 9 (MAP kinase 9, MAPK 9) JNK-55 Stress-activated protein kinase 1a (SAPK1a) Stress-activated protein kinase JNK2 c-Jun N-terminal kinase 2 MAPK9 JNK2 PRKM9 SAPK1A
Observed MW:	44-55 kDa
Uniprot ID :	P45984

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04448(JNK2 antibody) at dilution of 1:1000