

AIMP2 antibody

Product Information

Catalog No.:	FNab04453
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	JTV1 gene
Synonyms:	Aminoacyl tRNA synthase complex-interacting multifunctional protein 2 Multisynthase complex auxiliary component p38 Protein JTV-1 AIMP2 JTV1
Observed MW:	36 kDa
Uniprot ID :	Q13155

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IHC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab04453(JTV1 Antibody) at dilution of 1:50



IP Result of anti-JTV1 (IP:FNab04453, 3ug; Detection:FNab04453 1:1000) with HeLa cells lysate 800ug.



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04453(JTV1 antibody) at dilution of 1:1000