

KLK8 antibody

Product Information

Catalog No.:	FNab04460
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	kallikrein-related peptidase 8
Synonyms:	Kallikrein-8 (hK8) Neurosin (NP) Ovasin Serine protease 19 Serine protease TADG-14 Tumor-associated differentially expressed gene 14 protein KLK8 NRPN PRSS19 TADG14
Observed MW:	30 kDa
Uniprot ID :	O60259

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human ovary tumor tissue slide using FNab04460(Kallikrein 8 antibody) at dilution of 1:200



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04460(Kallikrein 8 antibody) at dilution of 1:500