

KATNAL1 antibody

Product Information

Catalog No.:	FNab04470
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Regulates microtubule dynamics in Sertoli cells, a process that is essential for spermiogenesis and male fertility. Severs microtubules in an ATP-dependent manner, promoting rapid reorganization of cellular microtubule arrays(By similarity).

Immunogen information


Immunogen:	katanin p60 subunit A-like 1
Synonyms:	Katanin p60 ATPase-containing subunit A-like 1 (Katanin p60 subunit A-like 1) p60 katanin-like 1 KATNAL1
Observed MW:	50-55 kDa
Uniprot ID :	Q9BW62

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-KATNAL1 (IP:FNab04470, 3ug;
Detection:FNab04470 1:500) with HeLa cells
lysate 3000ug.



HeLa cells were subjected to SDS PAGE followed
by western blot with FNab04470(KATNAL1
Antibody) at dilution of 1:300