

KAZALD1 antibody

Product Information

Catalog No.:	FNab04472
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in the proliferation of osteoblasts during bone formation and bone regeneration.
Promotes matrix assembly(By similarity).

Immunogen information

Immunogen:	Kazal-type serine peptidase inhibitor domain 1
Synonyms:	Kazal-type serine protease inhibitor domain-containing protein 1 KAZALD1
Observed MW:	None
Uniprot ID :	Q96I82

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab04472(KAZALD1 Antibody) at dilution of 1:200