

## KBTBD2 antibody

### Product Information

|               |  |
|---------------|--|
| Catalog No.:  | FNab04475  |
| Size:         | 100µg  |
| Form:         | liquid   |
| Purification: | Immunogen affinity purified  |
| Purity:       | ≥95% as determined by SDS-PAGE   |
| Host:         | Rabbit   |
| Clonality:    | polyclonal   |
| Clone ID:     | None   |
| IsoType:      | IgG  |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

### Background


None

### Immunogen information


|              |   |
|--------------|---|
| Immunogen:   | kelch repeat and BTB(POZ) domain containing 2   |
| Synonyms:    | Kelch repeat and BTB domain-containing protein 2 BTB and kelch domain-containing protein 1 KBTBD2 BKLHD1 KIAA1489 |
| Observed MW: | 60-70 kDa   |
| Uniprot ID : | Q8IY47  |

### Application

|                       |                                   |
|-----------------------|-----------------------------------|
| Reactivity:           | Human, Mouse, Rat                 |
| Tested Application:   | ELISA, IHC, WB                    |
| Recommended dilution: | WB: 1:500-1:2000; IHC: 1:20-1:200 |
| Image:                |                                   |



Immunohistochemistry of paraffin-embedded human colon cancer using FNab04475(KBTBD2 antibody) at dilution of 1:50



COLO 320 cells were subjected to SDS PAGE followed by western blot with FNab04475(KBTBD2 antibody) at dilution of 1:500