

## KCNB1-Specific antibody

### Product Information

Catalog No.:	FNab04480
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	potassium voltage-gated channel, Shab-related subfamily, member 1
Synonyms:	Potassium voltage-gated channel subfamily B member 1 Delayed rectifier potassium channel 1 (DRK1, h-DRK1) Voltage-gated potassium channel subunit Kv2.1 KCNB1
Observed MW:	~130 kDa
Uniprot ID :	Q14721

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded mouse brain tissue slide using FNab04480( KCNB1-Specific Antibody) at dilution of 1:200



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab04480( KCNB1-Specific Antibody) at dilution of 1:600