

KCNK5 antibody

Product Information

Catalog No.:	FNab04495
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


pH-dependent, voltage insensitive, outwardly rectifying potassium channel. Outward rectification is lost at high external K(+) concentrations.

Immunogen information

Immunogen:	potassium channel, subfamily K, member 5
Synonyms:	Potassium channel subfamily K member 5 Acid-sensitive potassium channel protein TASK-2 TWIK-related acid-sensitive K(+) channel 2 KCNK5 TASK2
Observed MW:	72 kDa
Uniprot ID :	O95279

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab04495(KCNK5 antibody) at dilution of 1:800