

KCNQ1DN antibody

Product Information

Catalog No.:	FNab04500
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	KCNQ1 downstream neighbor
Synonyms:	KCNQ1 downstream neighbor protein Beckwith-Wiedemann region transcript protein KCNQ1DN BWRT
Observed MW:	None
Uniprot ID :	Q9H478

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemical of paraffin-embedded human kidney using FNab04500(KCNQ1DN antibody) at dilution of 1:50