

## KCTD15 antibody

### Product Information

Catalog No.:	FNab04505
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)


### Background

#### Immunogen information


Immunogen:	potassium channel tetramerisation domain containing 15
Synonyms:	BTB/POZ domain-containing protein KCTD15 Potassium channel tetramerization domain-containing protein 15 KCTD15
Observed MW:	32 kDa
Uniprot ID :	Q96SI1

#### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:100
Image:	



Immunohistochemistry of paraffin-embedded mouse lung using FNab04505(KCTD15 antibody) at dilution of 1:100



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab04505(KCTD15 antibody) at dilution of 1:1000