

C2CD5 antibody

Product Information

Catalog No.:	FNab04535
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	KIAA0528
Synonyms:	C2 domain-containing protein 5 C2 domain-containing phosphoprotein of 138 kDa C2CD5 CDP138 KIAA0528
Observed MW:	115 kDa
Uniprot ID :	Q86YS7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:1000 ; IF: 1:10-1:100
Image:	



IP Result of anti-KIAA0528 (IP:FNab04535, 4ug;
Detection:FNab04535 1:700) with HeLa cells
lysate 520ug.



HeLa cells were subjected to SDS PAGE followed
by western blot with FNab04535(KIAA0528
antibody) at dilution of 1:1000