

MTUS2 antibody

Product Information

Catalog No.:	FNab04536
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	KIAA0774
Synonyms:	Microtubule-associated tumor suppressor candidate 2 Cardiac zipper protein Microtubule plus-end tracking protein TIP150 (Tracking protein of 150 kDa) MTUS2 CAZIP KIAA0774 TIP150
Observed MW:	150 kDa
Uniprot ID :	Q5JR59

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab04536(KIAA0774 Antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab04536(KIAA0774 Antibody) at dilution of 1:1000