

KHDC4 antibody

Product Information

Catalog No.:	FNab04538
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	KIAA0907
Synonyms:	KH homology domain-containing protein 4 Brings lots of money 7 Pre-mRNA splicing factor protein KHDC4 KHDC4 BLOM7 KIAA0907 SNORA80EHG
Observed MW:	68 kDa
Uniprot ID :	Q7Z7F0

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IHC, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using FNab04538(KIAA0907 Antibody) at dilution of 1:50



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab04538(KIAA0907 Antibody) at dilution of 1:300