

KIAA1217 antibody

Product Information

Catalog No.:	FNab04542
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Required for normal development of intervertebral disks.

Immunogen information

Immunogen:	KIAA1217
Synonyms:	Sickle tail protein homolog KIAA1217 SKT
Observed MW:	120 kDa
Uniprot ID :	Q5T5P2

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:1000
Image:	



mouse colon tissue were subjected to SDS PAGE followed by western blot with FNab04542(KIAA1217 Antibody) at dilution of 1:600