

## RELCH antibody

### Product Information

Catalog No.:	FNab04544
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	KIAA1468
Synonyms:	RAB11-binding protein RELCH LisH domain and HEAT repeat-containing protein KIAA1468 RAB11 binding and LisH domain, coiled-coil and HEAT repeat-containing RAB11-binding protein containing LisH, coiled-coil, and HEAT repeats RELCH KIAA1468
Observed MW:	100-135 kDa
Uniprot ID :	Q9P260

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab04544(KIAA1468 antibody) at dilution of 1:500