

## NCKAP5L antibody

### Product Information

Catalog No.:	FNab04547
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	KIAA1602
Synonyms:	Nck-associated protein 5-like (NCKAP5-like) Centrosomal protein of 169 kDa (Cep169) NCKAP5L CEP169 KIAA1602
Observed MW:	145 kDa
Uniprot ID :	Q9HCH0

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skin using FNab04547(KIAA1602 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04547(KIAA1602 antibody) at dilution of 1:300