

KIFC1 antibody

Product Information

Catalog No.:	FNab04573
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Minus end-directed microtubule-dependent motor required for bipolar spindle formation(PubMed:15843429). May contribute to movement of early endocytic vesicles(By similarity). Regulates cilium formation and structure(By similarity).

Immunogen information


Immunogen:	kinesin family member C1
Synonyms:	Kinesin-like protein KIFC1 Kinesin-like protein 2 Kinesin-related protein HSET KIFC1 HSET KNSL2
Observed MW:	55-70 kDa
Uniprot ID :	Q9BW19

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IF, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab04573(KIFC1 Antibody) at dilution of 1:50



rat testis tissue were subjected to SDS PAGE followed by western blot with FNab04573(KIFC1 Antibody) at dilution of 1:600