

KNG1 antibody

Product Information

Catalog No.:	FNab04577
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	kininogen 1
Synonyms:	Kininogen-1 Alpha-2-thiol proteinase inhibitor Fitzgerald factor High molecular weight kininogen (HMWK) Williams-Fitzgerald-Flaujeac factor Kininogen-1 heavy chain T-kinin Alternative names: Ile-Ser-Bradykinin Bradykinin Alternative names: Kallidin I Lysyl-bradykinin Alternative names: Kallidin II Kininogen-1 light chain Low molecular weight growth-promoting factor KNG1 BDK KNG
Observed MW:	60 kDa
Uniprot ID :	P01042

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human endometrial cancer using FNab04577(KNG1 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab04577(KNG1 antibody) at dilution of 1:300