

KIR3DX1 antibody

Product Information

Catalog No.:	FNab04584
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	killer cell immunoglobulin-like receptor, three domains, X1
Synonyms:	Putative killer cell immunoglobulin-like receptor-like protein KIR3DX1 Leukocyte receptor cluster member 12 KIR3DX1 KIR3DL0 LENG12
Observed MW:	None
Uniprot ID :	Q9H7L2

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab04584(KIR3DX1 Antibody) at dilution of 1:50