

KLF12 antibody

Product Information

Catalog No.:	FNab04598
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Confers strong transcriptional repression to the AP-2-alpha gene. Binds to a regulatory element(A32) in the AP-2-alpha gene promoter.

Immunogen information

Immunogen:	Kruppel-like factor 12
Synonyms:	Krueppel-like factor 12 Transcriptional repressor AP-2rep KLF12 AP2REP
Observed MW:	65 kDa
Uniprot ID :	Q9Y4X4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab04598(KLF12 Antibody) at dilution of 1:600