

KLHL42 antibody

Product Information

Catalog No.:	FNab04608
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	kelch domain containing 5
Synonyms:	Kelch-like protein 42 Cullin-3-binding protein 9 (Ctb9) Kelch domain-containing protein 5 KLHL42 KIAA1340 KLHDC5
Observed MW:	56 kDa
Uniprot ID :	Q9P2K6

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



mouse thymus tissue were subjected to SDS PAGE followed by western blot with FNab04608(KLHDC5 antibody) at dilution of 1:500