

KRAS-2B-specific antibody

Product Information

Catalog No.:	FNab04644
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	v-Ki-ras2 Kirsten rat sarcoma viral oncogene homolog
Synonyms:	GTPase KRas K-Ras 2 Ki-Ras c-K-ras c-Ki-ras GTPase KRas, N-terminally processed KRAS KRAS2 RASK2
Observed MW:	21 kDa
Uniprot ID :	P01116

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF, IHC, FC
Recommended dilution:	WB: 1:500 - 1:1000; IHC: 1:50 - 1:500; IF: 1:50 - 1:500
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab04644(KRAS-2B-specific antibody) at dilution of 1:50



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab04644(KRAS-2B-specific antibody) at dilution of 1:1000