

KRT12 antibody

Product Information

Catalog No.:	FNab04650
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a unique role in maintaining the normal corneal epithelial function. Together with KRT3, essential for the maintenance of corneal epithelium integrity(By similarity).

Immunogen information

Immunogen:	keratin 12
Synonyms:	Keratin, type I cytoskeletal 12 Cytokeratin-12 (CK-12) Keratin-12 (K12) KRT12
Observed MW:	50-55 kDa
Uniprot ID :	Q99456

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	



mouse eye tissue were subjected to SDS PAGE followed by western blot with FNab04650(KRT12 antibody) at dilution of 1:800