

KRT23 antibody

Product Information

Catalog No.:	FNab04651
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	keratin 23(histone deacetylase inducible)
Synonyms:	Keratin, type I cytoskeletal 23 Cytokeratin-23 (CK-23) Keratin-23 (K23) KRT23
Observed MW:	48 kDa
Uniprot ID :	Q9C075

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab04651(KRT23 Antibody) at dilution of 1:50



BxPC-3 cells were subjected to SDS PAGE followed by western blot with FNab04651(KRT23 antibody) at dilution of 1:1000