

KRT40 antibody

Product Information

Catalog No.:	FNab04655
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May play a role in late hair differentiation.

Immunogen information

Immunogen:	keratin 40
Synonyms:	Keratin, type I cytoskeletal 40 Cytokeratin-40 (CK-40) Keratin-40 (K40) Type I hair keratin Ka36 KRT40 KA36
Observed MW:	48 kDa
Uniprot ID :	Q6A162

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:200-1:1000
Image:	



A431 cells were subjected to SDS PAGE followed by western blot with FNab04655(KRT40 Antibody) at dilution of 1:500