

## KRT80 antibody

### Product Information

Catalog No.:	FNab04657
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	keratin 80
Synonyms:	Keratin, type II cytoskeletal 80 Cytokeratin-80 (CK-80) Keratin-80 (K80) Type-II keratin Kb20 KRT80 KB20
Observed MW:	47 kDa, 50 kDa
Uniprot ID :	Q6KB66

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab04657(KRT80 antibody) at dilution of 1:200



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with FNab04657(KRT80 antibody) at dilution of 1:1500