

KRT81 antibody

Product Information

Catalog No.:	FNab04658
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	keratin 81
Synonyms:	Keratin, type II cuticular Hb1 Hair keratin K2.9 Keratin, hair, basic, 1 Keratin-81 (K81) Metastatic lymph node 137 gene protein (MLN 137) Type II hair keratin Hb1 Type-II keratin Kb21 ghHKb1 (ghHb1) KRT81 KRTHB1 MLN137
Observed MW:	54 kDa
Uniprot ID :	Q14533

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNa04658(KRT81 antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNa04658(KRT81 antibody) at dilution of 1:300