

## LTF antibody

### Product Information

Catalog No.:	FNab04674
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	lactotransferrin
Synonyms:	Lactotransferrin (Lactoferrin) Growth-inhibiting protein 12 Talalactoferrin Lactoferricin-H (Lfcin-H) Kaliocin-1 Lactoferroxin-A Lactoferroxin-B Lactoferroxin-C LTF GIG12 LF
Observed MW:	51 kDa
Uniprot ID :	P02788

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab04674(LTF antibody) at dilution of 1:50



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab04674(LTF antibody) at dilution of 1:500